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Before using the microscopes, new users should attend the microscopy training workshop (Video) to ensure a full understanding of the basics of microscopy.

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Observation of operation by a certified user three times (a total of 3 h) is required (signatures of a new user are required).

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Operation under the supervision of a certified user for three times (a total of 3 h) is required (signatures of a certified user are required).

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After completing the aforementioned assignment, please email a signed Agreement Sheet of Instrument by the lab PI to the TA.

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Check all the records with the TA and take an operation exam. Then apply for online reservation password and ID-card approval for accessing room 927.

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Users who have accumulated 20 hours of experience and have no violation records may apply for nighttime and weekend access. Please contact Professor Li-Heng Tuan for authorization.

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To reserve time, please register on the microscope paper sign-in sheet and accurately record your usage hours. If a user reserves time but fails to show up three times, their access will be suspended for one month.

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Every current user should turn off the power supplies of all devices before leaving the room.

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If any user is found or reported to leave without turning off any of the devices, this will be considered misconduct.

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Please notify the administrator one week in advance if you wish to use the 100x objective lens. Do not replace the objective lens yourself.

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If you need to change the specimen holder, please revert to the standard slide holder after use.

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After use, turn the objective lens back to 4x and lower the stage to its lowest position.

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To ensure smooth operation, all users are prohibited from changing any setting in the microscope software.

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Please confirm the software remained in its default setting before use. If the user was found or reported to alter settings and not returning to default when finished, the last user will be held accountable and an offense will be registered to the user.

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Users are encouraged to save their settings in separate files or reuse the configuration of an existing picture.

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Users have to turn off the temperature control and CO₂ pump after usage and make sure the main valve of CO₂ tank is closing. The user who doesn't turn it off will record an offense.

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Users must book a time block each time before use. Users must also fill out the form in the notebook after using the microscope. Users have to cancel the reservation if users cannot use the microscope on time.

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Users without certification are not allowed to use the microscope.

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If the user only uses the computer for data processing without using the microscope, you must still record it in the log for reference. If you need to use the computer for more than 30 minutes, a reservation is required and there will be a fee.

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Pure image processing can use offline software on the neighboring computer.

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Please use a formatted external hard drive to transfer the data to avoid infecting the computer with viruses.

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Backing up data in a separate storage device immediately is strongly encouraged.

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Files will be removed from the hard drive periodically at the beginning of each month and will be announced in advance.

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If you find any abnormality in operating the microscope, you need to

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In case of emergency, call the administrator (Li-Heng Tuan) immediately.

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Users are responsible for cleaning the microscope room (including the floor, desks, the exterior of the microscopes, and oil immersion lenses) and taking out the trash.

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Do not eat or drink in the microscope room.

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Users will be subjected to one offense for disobeying the above rules.

III. 操作不當之處罰 Punishment for misconduct

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A verbal warning will be issued for the first offense.

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Individual users will be penalized for a one-month restriction of use for the second offense and the instructing professor will be notified.

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The whole lab will be penalized for a month for the third offense.

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The record of misconduct runs every year.

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Improper operation is defined as follows:

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Failure to turn off any equipment power before leaving will be counted as one improper operation.

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Unauthorized modification of any software default settings will be counted as one improper operation.

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For users of temperature control and CO₂ systems, failure to turn off the temperature control and CO₂ pump before shutting down will be counted as one improper operation.

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For users of CO₂, failure to close the main CO₂ cylinder valve tightly before leaving will be counted as one improper operation.

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Failure to clean oil-immersion objectives after use will be counted as one improper operation.

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Failure to return the objective lens to 4× after use will be counted as one improper operation.

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Failure to remove personal trash after use will be counted as one improper operation.

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Eating or drinking inside the microscope room will be counted as one improper operation.

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If an unqualified user uses any microscope alone, this user will be forbidden to use all fluorescence microscopes for one semester. Individual user's labs will be responsible for all collateral damages and held accountable for the offense.

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The maintenance fee and cost of consumable parts of the microscope (please see appendix) will be partitioned to each laboratory according to total use time.

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The bill will be emailed to both individual PI and secretary in April and October every year, and money will be deducted from PI's departmental account directly.

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Usage fee of THUNDER for PI in College of Life Science (see less than 30 min as 30 min)

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Usage fee of THUNDER for PI NOT in College of Life Science (see less than 30 min as 30 min)

Basic usage: NT 300/ 30 min

CO₂ tank: NT 75/ 30 min