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## Regulations and Rules for use of SMZ25 stereo microscope

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observe operation by a certified user three times, for a total of 3 hours (signatures of a new user are required)

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operate with a certified user aside three times, for a total of 3 hours (signatures of a certified user are required)

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check all the records with TA and take an operation exam. Then apply for online reservation password and ID-card approval for accessing room 332.

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Before using the microscopes, new users are encouraged to attend the microscopy training workshop or have taken any related course about microscopy to ensure a full understanding of the basics of microscopy,

## II. 操作不當之處罰 Punishment about misconduct

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A verbal warning will be issued for the first offense.

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Individual user will be penalized for a one month restriction of use for the second offense and instructing professor will be notified.

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The whole lab will be penalized for a month for the third offense.

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The record of misconducts runs in a yearly basis.

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Individual lab is responsible for the cost of any damages due to the abuse of the instrument.

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If an unqualified user uses any microscope alone, this user will be forbidden to use all fluorescence microscopes for one semester. Individual user's lab will be responsible for all collateral damages and held accountable for the offense.

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The current user should turn off everything before leaving the room.

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If any user is found or reported to leave without turning off lamp, this will be considered as one misconduct.

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To ensure smooth operation, all users are prohibited to change any setting in the microscope software.

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Users are encouraged to save their settings in separate files, or **reuse** the configuration of an existing picture.

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Please confirm the software remained in its default setting before use. If the user was found or reported to alter settings and not returning to default when finished, the last user will be held accountable and an offense will be registered to the user.

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Users must book time block each time before use. Users must also fill the form in the notebook after using the microscope. Users have to cancel the reservation if users cannot use the microscope on time.

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Users without certification are not allowed to use reservation system.

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Users will be subjected to one offense for disobeying the above rules.

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The users should also fill the form in the notebook even if using the computer only.

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1. 強烈建議各使用者在每次掃瞄後，將自己的實驗資料備份。  
Backing up data in a separate storage device immediately is strongly encouraged.
2. 管理者將於每月月初移除公用電腦內的資料。  
Files will be removed from the hard drive periodically at the beginning of each month, and will be announced in advance.

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If you find any abnormality in operating the microscope, you need to

1. 在登記表上註明問題。  
make a note in the notebook
2. 通知教學助理。  
report to TA
3. 若有緊急狀況，請立即通知管理者（黃建銘）。  
In case of emergency, call the administrator (Jian-Ming Huang) immediately.

#### **F. 清潔 cleaning**

1. 使用者需負責本顯微鏡室的清潔（包含地面、桌面、顯微鏡外部、油鏡鏡頭等）  
Users are responsible for the cleaning of the microscope room, including the floor, the desk, the external part of the microscope, and the oil objectives, etc..
2. 本顯微鏡室禁止飲食。  
Do not eat or drink in the microscope room.
3. 違反上述規範者，視為一次操作不當記錄。  
Users will be subjected to one offense for disobeying the above rules.

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1. 顯微鏡之維修與消耗性費用(請見附錄)將由使用實驗室按時間比例分攤。  
The maintenance fee and cost of consumable parts of microscope (please see appendix) will be partitioned to each laboratory according to total using time.
2. 每年4月及10月各收費一次。帳單會以電子郵件寄到各指導教授及院秘書，費用將直接由教授的系所經費中扣除。  
The bill will be emailed to both individual PI and secretary on April and October every year, and money will be deducted from PI's departmental account directly.

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Usage fee of SMZ25 for PI in college of life science (see less than 30 min as 30 min)

Basic usage	NT 50/ 30min
-------------	--------------